
PaX-i Plus™

PaX-i Plus™

User Manual Version 1.39

Model : PCH-30CS

- English
- Français
- Deutsch
- Italiano
- Español
- русский

For the European region

ENGLISH

Notice	i
Conventions in this Manual.....	ii
1. Imaging System Overview	1
1.1 Equipment Overview	2
1.2 Emergency Stop Switch	5
1.3 Exposure Switch.....	6
2. Getting Started	7
2.1 Turning on the Equipment	7
2.2 The Recommended PC Requirements	8
2.3 Running the Image Viewer (EzDent-i).....	9
2.4 Creating a New Patient Record	10
2.5 Initiating the Console Software	12
3. Acquiring PANO Images	13
3.1 Configuring Exposure Parameters.....	13
3.2 Patient Positioning.....	18
3.3 X-ray Exposure.....	34
3.4 Finishing the Scan	35
3.5 Checking the Captured Images	35
4. Acquiring CEPH Images (Optional).....	37
4.1 Configuring Exposure Parameters.....	37
4.2 Patient Positioning.....	41
4.3 X-ray Exposure.....	51
4.4 Finishing the Scan	52
4.5 Checking the Captured Images	52

Notice

Thank you for purchasing the **PaX-i Plus (Model: PCH-30CS)** extra-oral imaging system.

The user manual is a part of product.

This manual describes how to operate the **PCH-30CS** system. It is recommended that you thoroughly familiarize yourself with this manual to make the most effective use of this equipment.

Observe all cautions, safety messages and warnings which appear in this manual.

Due to constant technological improvement, the manual may not contain the most updated information and is subject to change without prior notice to the persons concerned. For further information not covered in this manual, please contact us at:

VATECH Co., Ltd.

Phone: +82-1588-9510

E-mail: gcs@vatech.co.kr

Website: www.vatech.com

Manual Name: PaX-i Plus (Model: PCH-30CS) User Manual

Version: 1.39

Publication Date: 2022-09

Copyright by © 2017 VATECH Co., Ltd.

All rights reserved.

The documentation, brand name and logo used in this manual are copyrighted.

No part of this manual may be reproduced, transmitted, or transcribed without the expressed written permission of the manufacturer.

We reserve the right to make any alterations which may be required due to technical improvement.

Conventions in this Manual

The following symbols are used throughout this manual. Make sure that you fully understand each symbol and follow the instructions accompanied.

To prevent physical injury and/or damage to the equipment, please observe all warnings and safety information included in this document.

	WARNING	Indicates information that should be followed with the utmost care. Failure to comply with a warning may result in severe damage to the equipment or physical injury to the operator and/or patient.
	CAUTION	Indicates a situation that demands prompt and careful action, a specific remedy, or emergency attention.
	IMPORTANT	Indicates a situation or action that could potentially cause problems to the equipment and/or its operation.
	NOTE	Emphasizes important information or provides useful tips and hints.
	RADIATION	Indicates a possible danger from exposure to radiation.
	SINGLE USE	Indicates a component which must be replaced for each new patient.
	ESD susceptibility	Indicates that an item is susceptible to damage from electrostatic discharges.

1. Imaging System Overview

PaX-i Plus (Model: PCH-30CS) is an advanced 2-in-1 digital X-ray imaging system that incorporates PANO and CEPH (Optional) imaging capabilities into a single system.

The **PaX-i Plus** is intended for use to take a panoramic and cephalometric image of the oral and craniofacial anatomy to provide diagnostic information for adult and pediatric patients.

PaX-i Plus may only be used by dentists, X-ray technicians and other professionals who are licensed to perform X- rays by the law of the region in which it is being used.

Standards and Regulations

PaX-i Plus (Model: PCH-30CS) is designed and developed to comply with the following international standards and regulations:

- MEDICAL - APPLIED ELECTROMAGNETIC RADIATION EQUIPMENT
AS TO ELECTRICAL SHOCK, FIRE, AND MECHANICAL
HAZARDS ONLY IN ACCORDANCE WITH
ANSI/AAMI ES60601-1 (2005) + AMD 1 (2012), CAN/CSA-C22.2 No. 60601-1 (2014), IEC 60601-1-3 (2008), IEC 60601-2-63 (2012)
- 21 CFR 1020.30, 31, 33
- NEMA Standard publication PS 3.1-3.18


This is Class IIb equipment and obtained CE marking in April 2007 for regulations compliance in accordance with the revised European Union's MDD (Medical Devices Directive) 93/42 EEC.

Classifications (IEC 60601-1 6.1)

- Degree of protection against water ingress: Ordinary Equipment: IPX0
- Degree of protection against electric shock: Class 1 equipment, Type B Applied Parts (Chinrest, Bite and cover, Nasal Positioner and cover, Ear Rod and cap, Carpus Plate)


1.1 Equipment Overview

No.	Item	Description
1	X-ray Detector for CEPH (Optional)	Xmaru2602CF
2	Nasal Positioner	- Positions the patient during CEPH imaging. - The ruler used for reference in an acquired image that is different from actual size
3	Ear Rods	Secure the patient's head during CEPH imaging.
4	Emergency Stop Switch	Immediately stops the moving parts and cuts off all power to the equipment's electrical components.
5	D-Sub Connector	The input signal port for Column UP/DOWN Switch
6	Main Power Switch	Turns on / off the main power of the equipment.
7	Horizontal Beam Lever	Aligns the Horizontal Beam in PANO mode.
8	Handle Frame	Held firmly by the patient during imaging to stabilize his / her position.

No.	Item	Description
9	Column UP/DOWN Switch (optional)	Adjusts the height of the Vertical Frame.
10	Telescopic Column	Moves by the Column UP/DOWN button or switch for patient positioning.
11	Stationary Column	Supports whole part of the equipment.
12	Base (Optional)	Balances the equipment and maintains its safety.
13	Vertical Frame	Holds the Rotating Unit. Can be controlled with Column UP/DOWN switch.
14	X-ray Generator	The vacuum tube where the X-ray is produced.
15	Rotating Unit	Rotates around the patient's head while image is being acquired. (Its movement is different according to the scan mode.)
16	LED Lamp	Displays the status of X-ray exposure. - Green: Standby - Yellow: In operation
17	X-ray Detector for PANO	Xmaru1501CF-PLUS
18	Temple Supports	Supports patient's head by holding the temples. Used in the PANO mode.
19	Chinrest	The place to rest the chin.
20	Canine Teeth Beam Lever	Aligns the Canine Teeth Beam in PANO mode.
21	Temple Supports OPEN/CLOSE Wheel	Adjusts the Temple Supports for patient positioning.
22	Control Panel	Operates the Laser Beam and adjusts the height of the Vertical Frame. (For the details, refer to 1.1.1 Control Panel .) 

1.1.1 Control Panel

No.	Item	Description
1	Column UP/DOWN button	Moves the Vertical Frame up or down. (For adjusting the height of the Chinrest)
2	Laser Beam ON/OFF button	Turns the Laser Beam on / off.

1.2 Emergency Stop Switch

During operation, the following emergency situations may occur:

- X-ray emission even after the **Exposure Switch** has been released
- Physical injury to the patient or damage to the equipment
- Other emergency situations

If a problem occurs during image acquisition, press the red **Emergency Stop Switch** to immediately stop the moving parts and cut off all power to the equipment's electrical components. To restart the equipment, turn the **Emergency Stop Switch** clockwise until it pops up.

The **Emergency Stop Switch** is located under the bottom of the Vertical Frame.


1.3 Exposure Switch

The **Exposure Switch** allows the operator to control image acquisition from outside of the X-ray room.

Press and hold the **Exposure Switch** down until acquisition is completed. Premature release of the **Exposure Switch** will abort image acquisition.

Pressing the **Exposure Switch** activates the LED indicator to turn yellow. This color indicates that the X-ray is being emitted.

IMPORTANT

- The **Exposure Switch** is detachable. Ensure the **Exposure Switch** cable is not detached from the unit accidentally during the operation.
- Keep vocal/visual contact with the patient during exposure. If any problem occurs during exposure, release the **Exposure Switch** immediately.

2. Getting Started

2.1 Turning on the Equipment

WARNING

- Do not place the patient close to the equipment when it's being turned on. Doing so may cause physical injury to the patient and damage to the equipment.
- Do not operate the PC while the equipment is in operation. Doing so may cause an error on the equipment.

CAUTION

- Extreme fluctuation of temperature may cause condensation inside the equipment. Do not switch on the equipment until it has reached normal room temperature.
- Rebooting the equipment: After turning it off, wait for approx. 20 seconds before turning it on again.
- Warm-up the equipment for at least 5 minutes before the operation. For the best image quality, it is recommended to have a warm-up phase for more than 30 minutes.

IMPORTANT

If the equipment has not been used for a long time, please let it have enough time to be warmed up. It extends the life of the X-ray tube.

The imaging system mainly consists of the imaging equipment and the PC.

Before turning on the equipment, please confirm that the equipment and PC have been installed correctly.

- Turn on the PC.
- Press the **Main Power Switch** that is located under the bottom of the Vertical Frame to turn on the equipment.

NOTICE

Main Power Switch isolates its circuits electrically from the supply mains on all poles simultaneously.

- Make sure that the green LED light at the top of the equipment is on.

2.2 The Recommended PC Requirements

IMPORTANT

- It is mandatory to ensure that the PC system configuration is compatible with the PC system requirements for the imaging and image viewer software.
- Since image quality may be deteriorated from lack of resources, observe the requirement guideline specified in the following tables.
- The PC components shall be approved by UL/CSA.
- The PC shall be grounded well protectively.
- The multiple portable socket-outlets shall not be placed on the floor.
- In case the equipment is to be installed in area with an unstable electric power supply, it is strongly recommended to use the AVR (automatic voltage regulator) to keep the line voltage stable.

Item	Specifications
CPU	Intel® Core™ i3-7100 3.9GHz 3MB cache, 2cores
Chipset	Intel® Q270 Chipset
RAM	2X4GB DDR4-2400 DIMM NECC UNB
Hard disk drive	500GB SATA 7200 1st HDD
Graphic board	Integrated Intel® HD 630 Graphics
Ethernet interface	Integrated Intel® I219LM Gigabit Network Intel Ethernet I210-T1 PCIe x1 Gb NIC(Optional)
Serial Port (RS232)	1 (On board)
Power supply	≥ 400 Watts (93% Efficiency)
Slots	2 PCI Express x 1 Slot 2 PCI Express x 16 Slots
CD/DVD drive	DVDRW
Operating system	Windows 10 Professional 64 bit

IMPORTANT

- In Windows 10, disable Windows Defender.
-  When Windows Defender is not enabled, Windows 10 is not protected from malware and virus.

WARNING

- If you need to install a grabber card, plug it into the x4 express slot.

2.3 Running the Image Viewer (EzDent-i)

The Imaging Program is interfaced with **EzDent-i** and the user can analyze the image acquired from the Console Software easily and rapidly. On your desktop, double-click **EzDent-i** icon. The **EzDent-i** main window will be displayed.

NOTICE

For further details on this subject, refer to the **EzDent-i User Manual**.

NOTICE

Security Capabilities


- It is recommended to install and operate **EzDent-i** SW within a secure operating environment that allows only authorized users to access and a system network equipped with Windows built-in firewall, Windows Defender antispyware tools and other commonly used 3rd party security tools and application systems.
- The latest updates for anti-virus software and a firewall are recommended.
- The software can be updated by the manufacturer only. Unauthorized software update through a third party, not the manufacturer, is strictly prohibited. For cyber security issues related to the software and medical devices, please contact the manufacturer.

NOTICE

For **PCH-30CS** dental X-ray system, the Console Software is being accessed through 2D viewer (**EzDent-i**) SW. The Console Software does not have image storage capability of their own and will not be able to keep patient information.

2.3.1 Creating a New Patient Record

To create a new patient record with **EzDent-i**:


The screenshot shows the EzDent-i Main Window with the **PATIENT** tab selected. The window features a search bar, a recent patients list, and a table of patient records. The table lists patient information including Chart No., Name, and Date of Birth.

Chart No.	Name	Date of Birth
20191029_103203	Koo HyunJun	1974-08-15
20191029_103621	HyunJun Koo	1974-08-15
20191029_103722	Hyun koo	2019-01-01
202191029_103203	Koo HyunJun	1974-08-15

EzDent-i Main Window

1. Select the **PATIENT** tab and click . After clicking the button, a dialog box appears as the below.

The screenshot shows the ADD PATIENT dialog box. The dialog box contains fields for Chart No., Name, Gender, Date of Birth, Doctor 1, Doctor 2, Social ID, Phone, Mobile, E-Mail, Zip Code, and Address. The 'Add' button is highlighted.

2. In the dialog box, fill out **Chart No.** and **Name**. If necessary, fill out other areas such as social ID and mobile.
3. Click  to save the new patient record.


NOTICE

The chart number and name are required fields. You cannot leave them as blanks.

2.3.2 Retrieving Patient Records

To retrieve the patient record saved in **EzDent-i**:


1. Go to the search box and enter the **patient's name** or **chart number**.


NOTICE

You can use a virtual keyboard to search for the patient record. Click the keyboard icon next to the search box to use the virtual keyboard:

2. Click .
3. Check and double-click the patient record with the matching chart number or name.


2.3.3 Initiating the Console Software

To initiate the console software:


NOTICE

Before starting the console software, you must create new patient information. Go to **2.3.1 Creating a New Patient Record** for more information.

1. Search and select the patient record.
2. Click the **ACQUISITION** tab.

3. Select a modality you want to capture the image with. The displayed options can be different depending on your equipment's option.


3. Acquiring PANO Images

3.1 Configuring Exposure Parameters

To acquire PANO Images, **2. Getting Started** must be completed first. If not, you must return to the **2. Getting Started** and finish the step first.


NOTICE

You can set the imaging parameters on the Console Software running on the PC.

3. Acquiring PANO Images


1. Click PANO button in the Modality Option of the Main Screen.

NOTICE


“CEPH” button exists only when each imaging program is included in the equipment.

2. Select a Pano Option.

Mode	Description
Normal (Default)	- Provides a normal panoramic image.
Auto Focusing (Optional)	- Selectively provides specific multiple panoramic images having different focal planes.

3. Select an Image Option.

Mode	Description
UHD (Optional)	Ultra High Definition image
HD (Default)	High Definition image
Normal	Normal quality image

4. Make an Arch Selection.

Arch Selection	Description
Narrow	Panoramic image of V-shaped palatal arches (small number of adult females)
Normal	Panoramic image of normal adult palatal arches
Wide	Panoramic image of square-shaped palatal arches (some number of adult males)
Child	Panoramic image of child palatal arches, approximately more than 40% less x-ray dose than in Normal mode.
Orthogonal	Panoramic image where the x-ray angle enters vertically in between the teeth so overlapping images are minimized. NOTICE If Orthogonal Arch is selected, Bitewing examinations (Bitewing, Bitewing Incisor (Optional), Bitewing Right, Bitewing Left) are activated.

5. Select an Examination Program in the Pano Examination or Special Examination panel.

NOTICE

- To activate Bitewing examination options - Bitewing, Bitewing Incisor (Optional), Bitewing Right, Bitewing Left, select Orthogonal Arch in the Arch Selection panel.

- When a Special Examination option is clicked, "PANO Examination" panel is disabled. If you want to select a PANO Examination option, please conduct Arch selection again.

- The Gender / Age group of the patient is selected automatically based on the patient information. If necessary, you can select the option manually.

NOTICE

Gender / Age Group		VATECH's Standard
Child		2 ~ 12 years of age
Adult	Man	> 12 years of age
	Woman	

- Select X-ray intensity.

NOTICE

Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft :

Soft ≤ Normal ≤ Hard

Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Intensity
Child	53±3	>53±3	Hard
		53±3	Normal
		<53±3	Soft
Adult	56±3	>56±3	Hard
		56±3	Normal
		<56±3	Soft

8. The values of tube voltage and current are configured automatically according to the patient's gender / age group and X-ray intensity. Click the **UP/DOWN** arrow to adjust kVp and mA. The dose is adjustable by ± 1.0 kVp and ± 1.0 mA respectively.

9. Click **CONFIRM** button when exposure parameter setting is completed.

CONFIRM

NOTICE

When you click **CONFIRM** button,

- The Rotating Unit will move to its scanning position.
- **READY** button will be activated. (This means that the equipment is ready for X-ray exposure.)
- Three Laser Beams (Vertical Beam, Horizontal Beam and Canine Teeth Beam) will be activated.
 - The Laser Beams are turned off automatically after 20 minutes or when **READY** button is clicked.
- The DAP (Dose Area Product), Scan Time and Exposure Time will be displayed below Patient Information window.

DAP
127.334307 mGy x cm²

Scan-time
13.5 Sec

Exposure-time
13.5 Sec

10. Guide the patient to the equipment.

3.2 Patient Positioning


- Have patients (especially pregnant women and children) wear a lead apron to protect themselves from residual radiation.
- Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.

- Correct posture reduces the shadow casted by the patient's cervical spine and allows clear image acquisition.
- Metal implants or bridges may reduce the quality of the images.
- Be sure to adjust laser beam correctly. Otherwise, the quality of images can be lower due to ghost images or expansion / reduction of the images.

Getting prepared

1. Let the patient remove all the metal objects (glasses, earrings, hair pins, braces, false teeth, etc.). Metal objects may induce ghost images and lower the image quality.
2. Have the patient wear a lead apron to protect themselves from residual radiation.
3. Use the **Column UP/DOWN** button or switch option to adjust the equipment to match the height of the patient.
4. Press the Rotation test button in the Console program and check that the equipment touches a part of the patient while the equipment is running.


- In general, imaging is performed with the patient in an upright position. However, a stool may be used for imaging patients with special circumstances. If a stool is used, ensure that the beams and movement of unit are not obstructed by the stool.

3.2.1 PANO Examination Mode (Standard / Right / Left / Front / Orthogonal)

Normal Patient Positioning

1. Insert the Normal Bite into the Normal Chinrest and cover it with a Sanitary Vinyl Cover.


- The Sanitary Vinyl Cover is for single use only. It should be replaced for each patient. Be sure to use the approved vinyl cover.

- Clean the Chinrest and the Bite with ethanol and wipe with a dry towel before the next patient.


2. Loosen the **Temple Supports OPEN/CLOSE Wheel** under the control panel to widen the Temple Supports.

Temple Supports
OPEN/CLOSE
Wheel


3. Guide the patient to the inside of the equipment.
4. Use the **Column UP/DOWN** button or switch option to adjust the height of the equipment so that the patient's chin reaches the Chinrest.

Column
UP/DOWN
button


Column
UP/DOWN
switch



3. Acquiring PANO Images

5. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
 - Hold the handles tightly.
 - Press the chest against the equipment.
 - Keep both feet close to the inside of the base.
 - Keep both shoulders parallel.
 - Straighten the Cervical Spine and stand still.
6. Let the patient bite the Bite along its grooves with his/her front teeth.

7. Let the patient maintain the posture as follows:
 - Close the mouth.
 - Place the tongue to the roof of the mouth.
 - Close the eyes.


To acquire the best image possible, ask the patient not to:

- Breathe or swallow saliva during image acquisition
- Move during image acquisition


Edentulous Patient Positioning

1. Insert the **Edentulous Bite** into the Normal **Chinrest**.


Clean the Chinrest and the Bite with ethanol and wipe with a dry towel before the next patient.


2. Loosen the **Temple Supports OPEN/CLOSE Wheel** under the control panel to widen the Temple Supports.

3. Guide the patient to the equipment.
4. Use the **Column UP/DOWN** button or switch option to adjust the height of the equipment so that the patient's chin reaches the Chinrest.

Column
UP/DOWN
switch

5. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
 - Hold the handles tightly.
 - Press the chest against the equipment.
 - Keep both feet close to the inside of the base.
 - Keep both shoulders parallel.
 - Straighten the Cervical Spine and stand still.
6. Let the patient maintain the posture as follows:
 - Close the mouth.
 - Place the tongue to the roof of the mouth.
 - Close the eyes.

Laser Beam Aligning

Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.



If the laser beam is not correctly positioned, there may be distortion, causing the image to be enlarged or reduced, or ghost shadows may occur and lower the image quality. Be sure to align the laser beam properly.

- Three Laser Beams (Vertical Beam, Horizontal Beam and Canine Teeth Beam) will be activated when **CONFIRM** button is clicked.
- The Laser Beams are turned off automatically after 20 minutes or when **READY** button is clicked.)


- To turn the Laser beams on/off manually, click the  icon on the Control Panel of the Handle Frame or the  icon on the Console Software.

- The Vertical Beam is fixed. Align the center of the patient's face (Mid-sagittal Line) with the Vertical Beam. (It's to prevent horizontal expansion of the image)
- Align the Horizontal Beam in a straight line to the Frankfurt Line on the patient's face. Use the **Horizontal Beam** lever on the column (left side of the Control Panel) to position it. Make sure that the Horizontal Beam is aligned to the patient's face horizontally.


3. Acquiring PANO Images

3. Direct the patient to smile and align the Canine Teeth Beam to the center of the canines. Use the Canine Teeth Beam Lever (left side of the Control Panel) to adjust the position of the beam.

Finishing Patient Positioning

1. After checking the positions of the patient and the laser beam, tighten the **Temple Supports OPEN/CLOSE Wheel** under the control panel to prevent the patient's head from moving.

2. Click the **READY** button on the Console Software. The X-ray exposure has not started yet.

Make sure that the Temple Supports are in closed position before clicking the **READY** button.

3. Now go to **3.3 X-ray Exposure** to start the exposure.

3.2.2 SPECIAL Examination Mode (TMJ / Sinus)

3.2.2.1 TMJ Open Mode (LAT / PA)

You need to take three steps to acquire a TMJ open image: 1) positioning patient, 2) aligning laser beam and 3) X-ray exposure. You must complete the procedure for the TMJ open mode first before taking a TMJ close image.

Follow the steps below to acquire a TMJ open image.


Step 1: Patient Positioning

1. Remove the **Normal Chinrest** and Insert the **SINUS/TMJ bite** into the chinrest receptacle.

Ensure to clean Sinus/TMJ bite with ethanol and dry it with a towel each time before taking an image of patient.


2. Turn the **Temple Support Wheel** under the control panel to open the Temple Support wide for the patient.

3. Guide the patient to the inside of the equipment.

3. Acquiring PANO Images

4. Press the **Column UP/DOWN** button on the control panel or switch to adjust the height of the equipment until you see the patient's chin touches on the SINUS/ TMJ bite.


5. Ask the patient to stand in the center of the equipment and to do the following:
 - Grab the handles on both sides
 - Press the chest against the equipment
 - Place feet inside of the base.
 - Keep both shoulders parallel.
 - Stand upright and stretch the neck straight
6. Ask the patient to press the acanthion against the **SINUS/TMJ bite** and tilt the head forward by 5°.

IMPORTANT

Ensure that the patients do not touch the equipment with their jaws to maintain the proper position.

7. Ask the patient to maintain the posture as below, until the scanning is completed :
 - Keep the mouth open
 - Keep both eyes closed
 - Place the tongue on the roof of the mouth
 - Breath with the nose

Step 2: Laser Beam Alignment

DO NOT shine the laser beam directly into the person's eyes. This can cause vision loss or other serious damages to the eyes.



Ensure that patient's position is properly aligned with laser beams before starting an X-ray exposure. Wrong positioning of a laser beam or person can create shadows or distortions in the image.

- When you click the **CONFIRM** button, three laser beams (vertical, horizontal, and canine teeth) are activated at once.
- All beams are turned off automatically after 20 minutes or when you click the **READY** button.
- If you want to turn on or off the laser beams manually, click the  icon on the Control Panel of the Handle Frame or the  icon on the Console Software.


1. The Vertical Beam is fixed. Align the center of the patient's face (Mid-sagittal Line) with the Vertical Beam. (It's to prevent the horizontal expansion of the image)
2. Align the Horizontal Beam in a straight line to the Frankfurt Line on the patient's face. Use the **Horizontal Beam** lever on the column (left side of the Control Panel) to position it. Make sure that the Horizontal Beam is aligned to the patient's face horizontally.

Step 3: Finishing Patient Positioning

1. After checking the positions of the patient and the laser beam, turn the Temple Support Wheel again to close temple supports to keep the patient's head from moving.

Temple Support Wheel

2. Click the **READY** button on the Console Software. X-ray exposure has not started yet.

Ensure that the Temple Supports are in the closed position before clicking the **READY** button.


3. Now go to **3.3 X-ray Exposure** to start the exposure.

3.2.2.2 TMJ Close Mode (LAT / PA)

When you completed the procedure for TMJ open mode, the system will be ready for the TMJ close mode. Follow the steps below to acquire a TMJ close image.

Step 1: Patient Positioning

1. When you see the message "Do you want to capture a TMJ Close image?" Click the **OK** to start the TMJ close mode.


2. Turn the **Temple Support Wheel** under the control panel to open the temple support wide for the patient.
3. Guide the patient to the inside of the equipment.
4. Ask the patient to stand in the center of the equipment and to do the following:
 - Grab the handles on both sides
 - Press the chest against the equipment
 - Place feet inside of the base.
 - Keep both shoulders parallel.
 - Stand upright and stretch the neck straight
5. Ask the patient to press the acanthion against the **SINUS/TMJ bite** and tilt the head forward by 5°.

IMPORTANT

Ensure that the patient does not touch the equipment with their jaws to keep the proper position.

6. Ask the patient to maintain the posture as follows until the scanning is completed:
 - Close the mouth
 - Place the tongue on the roof of the mouth
 - Close eyes
 - Breathe with nose

NOTICE

The support unit of the integrated Chinrest must touch the patient's acanthion.

Step 2: Laser Beam Aligning

This is the same as the one for TMJ Open mode.

Step 3: Finishing Patient Positioning

This is the same as the one for TMJ Open mode.

3.2.2.3 Sinus Mode (LAT / PA)

Follow the steps below to acquire an image on Sinus mode.


Step 1: Patient Positioning

1. Remove the **Normal Chinrest** and Insert the **SINUS/TMJ bite**.


Clean the Chinrest and the Bite with ethanol and wipe with a dry towel before the next patient.


2. Turn the **Temple Supports Wheel** under the control panel to open the Temple Supports wide for the patient.

3. Guide the patient to the inside of the equipment.
4. Press the **Column UP/DOWN** button on the control panel or switch to adjust the height of the equipment until you see the patient's chin touches on the **SINUS/TMJ bite**.

5. Guide the patient to stand in the center of the equipment and direct them to remain in the position outlined below.
 - Grab the handles on both sides
 - Press the chest against the equipment
 - Place feet inside of the base
 - Keep both shoulders parallel
 - Stand upright and stretch the neck straight
6. Guide the patient to press the acanthion against the SINUS/TMJ bite and tilt the head backward about $10^{\circ} \sim 15^{\circ}$.
7. Ask the patient to maintain the posture as below until the scanning is completed :
 - Close the mouth
 - Place the tongue on the roof of the mouth
 - Close eyes
 - Breath with nose and stop swallowing

Step 2: Laser Beam Alignment**WARNING**

Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.



**CAUTION**

If the laser beam is not correctly positioned, there may be distortion, causing the image to be enlarged or reduced, or ghost shadows may occur and lower the image quality. Be sure to align the Laser Beam properly.

1. The Vertical Beam is fixed. Align the center of the patient's face (Mid-sagittal Line) with the Vertical Beam. (It's to prevent horizontal expansion of the image)
2. Align the Horizontal Beam in a straight line to the Frankfurt Line on the patient's face. Use the **Horizontal Beam** lever on the column (left side of the Control Panel) to position it. Make sure that the Horizontal Beam is aligned to the patient's face horizontally.

Vertical Beam
/ Mid-sagittal Line

Horizontal Beam
/ Frankfurt Line


Horizontal Beam Lever

Step 3: Finishing Patient Positioning

1. After confirming that both patients and laser beams are properly positioned, turn the **Temple Support Wheel** to close the part to keep the patient's head from moving.

Temple Support Wheel

2. Click the **READY** button on the Console Software. X-ray exposure has not started yet.

**CAUTION**

Ensure that the Temple Supports are in the closed position before clicking the **READY** button.

3. Now go to **3.3 X-ray Exposure** to start the exposure.

3.3 X-ray Exposure

- If an emergency occurs during image acquisition, release the **Exposure Switch** to cease X-ray emission.
- The operator shall observe the X-ray safety regulations applicable to his/her area at all times during the operation of this equipment.


- The operator must keep vocal/visual contact with the patient at all times during image acquisition process.
- Do not operate the PC during exposure. Doing so may cause the system to malfunction.

- Let the patient close the eyes during the operation.
- To acquire optimized images, instruct the patient to hold his/her breath and not to swallow. Also don't let the patient move until the Temple Supports are open.

1. Get out of the X-ray room and close the door.
2. Press and hold down the **Exposure Switch** until image acquisition is completed.

Yellow
- In Operation (X-ray On)

The image appears on the screen.

During X-ray exposure, the status appears as follows.

- The LED light of the **Exposure Switch** turns yellow.
- The LED light on the top of the equipment turns yellow.
- An alert sound comes out to indicate that X-ray emission is currently underway.
- On Console Software, the radiation mark turns yellow and "X-RAY" changes to "X-RAY ON".

X-RAY ON

3. Release the Exposure Switch when "Image capturing is completed" message appears on the screen.

3.4 Finishing the Scan

1. Open the Temple Supports and guide the patient out of the equipment.
2. For Normal Bite, remove the Sanitary Vinyl Cover from the Bite.
3. Press **READY** button on the Console Software to bring the Rotating Unit back to its initial position.

3.5 Checking the Captured Images

Acquired images can be reconstructed and converted to DICOM format.

The exported images can be confirmed in **EzDent-i**.

NOTICE

Refer to the **EzDent-i** User Manual for more information.

1. The images are transferred to **EzDent-i** automatically.
2. The images are automatically saved if automatic save option is configured as default. If it is not configured as default, click **Save** button to save the images.
3. To check the image, double-click the one on the **Patient List**.

This Page Left Blank Intentionally


4. Acquiring CEPH Images (Optional)

4.1 Configuring Exposure Parameters

To acquire CEPH images, **2. Getting Started** must be completed first.


NOTICE

You can set the imaging parameters on the Console Software running on the PC.

4. Acquiring CEPH Images (Optional)


1. Click **CEPH** button in the Main Screen.

NOTICE

“**CEPH**” button exists only when each imaging program is included in the equipment.

2. Select an Image Option.


Mode	Description
Normal	Normal quality image
Fast	Low dose image

NOTICE

Depending on the patient information, the **Image Option** is automatically selected as **Normal** for adult(Man / Woman) and **Fast** for child.

3. Select an examination program in the CEPH Examination panel.

4. The Gender / Age group of the patient is selected automatically based on the patient information. If necessary, you can select the option manually.

NOTICE

Gender / Age Group		VATECH's Standard
Child		2 ~ 12 years of age
Adult	Man	> 12 years of age
	Woman	

5. Select X-ray intensity.

NOTICE

Depending on the circumference of the patient's head, X-ray intensity may be classified as Hard, Normal, or Soft :

Soft ≤ Normal ≤ Hard

Age Group	Average Head Circumference (cm)	Range (cm)	X-ray Intensity
Child	53±3	>53±3	Hard
		53±3	Normal
		<53±3	Soft
Adult	56±3	>56±3	Hard
		56±3	Normal
		<56±3	Soft

6. The values of tube voltage and current are configured automatically according to the patient's gender / age group and X-ray intensity. Click the **UP/DOWN** arrow to adjust kVp and mA. The dose is adjustable by ±1.0 kVp and ±1.0 mA respectively.

- Click **CONFIRM** button when exposure parameter setting is completed.

CONFIRM

NOTICE

When you click **CONFIRM** button,

- **READY** button will be activated. (This means that the equipment is ready for X-ray exposure.)
- The DAP (Dose Area Product), Scan Time and Exposure Time will be displayed below Patient Information window.

DAP
127.334307 mGy x cm²

Scan-time
13.5 Sec

Exposure-time
13.5 Sec


- Guide the patient to the equipment.

4.2 Patient Positioning

 WARNING	- Have patients (especially pregnant women and children) wear a lead apron to protect themselves from residual radiation. - Be careful not to shine the laser beam directly into the person's eyes. Doing so may result in vision loss.
 CAUTION	Ensure that the Nasal Positioner is left unfolded, before adjusting the Ear Rods in the proper direction.
IMPORTANT	- Correct posture reduces the shadow casted by the patient's cervical spine and allows clear image acquisition. - Metal implants or bridges may reduce the quality of the images.

Getting prepared

1. Let the patient remove all the metal objects (glasses, earrings, hair pins, braces, false teeth, etc.). Metal objects may induce ghost images and lower the image quality.
2. Have the patient wear a lead apron to protect themselves from residual radiation.
3. Use the **Column UP/DOWN** button or switch option to adjust the equipment to match the height of the patient.

NOTICE

In general, imaging is performed with the patient in an upright position. However, a stool may be used for imaging patients with special circumstances. If a stool is used, ensure that the beams and movement of unit are not obstructed by the stool.

4.2.1 Lateral / Full Lateral Mode

NOTICE


Correct posture reduces the shadow casted by the patient's cervical spine and allows clear image acquisition.

Patient Positioning

1. Turn the Nasal Positioner to the **Lateral** mode Positioning Marker as shown below.

2. Leave enough space between the Ear Rods.

NOTICE


Use the Porion Position Reference Indicator that appears in the acquired image to easily confirm the location of Porion.

3. Guide the patient to the CEPH unit.
4. Direct the patient to relax his/her neck and shoulders and stand upright.
5. Use the **Column UP/DOWN** button or switch option to adjust the height of the CEPH unit to approximately match the height of the patient.


WARNING

After adjusting the height of the column, align the Ear Rods and Nasal Positioner to the patient.

6. Align the Ear Rods to the patient's ears properly so that the head does not move during the operation. And align the Nasal Positioner to the patient's nasion by adjusting its height.

7. Align horizontally so the patient's Frankfurt Line is parallel with the floor.

8. Direct the patient to swallow first before closing the mouth and to remain in his/her current position until image acquisition is completed.
9. Click the **READY** button on Console Software. No X-ray will be emitted at this point.
10. Now go to **4.3 X-ray Exposure** to start the exposure.


4.2.2 PA Mode

Patient Positioning

1. Turn the Nasal Positioner to the **PA / Waters' view / Carpus** mode Positioning Marker as shown below.

2. Fold the Nasal Positioner up. The Nasal Positioner is not used in PA mode.
3. Guide the patient to the CEPH unit.
4. Ask the patient to stand upright facing the X-ray detector. Make sure that the patient's shoulders are level and that his/her neck is relaxed.


5. Use the **Column UP/DOWN** button or switch option to adjust the height of the CEPH unit to approximately match the height of the patient.

WARNING After adjusting the height of the column, align the Ear Rods to the patient.

6. During the operation, properly align the Ear Rods to the patient's ears so his/her head does not move.
7. Align horizontally so the patient's Frankfurt Line is parallel with the floor.

8. Direct the patient to swallow first before closing his/her mouth and to remain in his/her current position until image acquisition is completed.
9. Click the **READY** button on Console Software. No X-ray will be emitted at this point.
10. Now go to **4.3 X-ray Exposure** to start the exposure.


4.2.3 SMV Mode

Patient Positioning

1. Turn the Nasal Positioner to the **SMV** mode Positioning Marker as shown below.

2. Fold the Nasal Positioner up. The Nasal Positioner is not used in SMV mode.
3. Guide the patient to the CEPH unit.
4. Guide the patient to face the X-ray tube and stand upright.


5. Use the **Column UP/DOWN** button or switch option to adjust the height of the CEPH unit to approximately match the height of the patient.

After adjusting the height of the column, align the Ear Rods to the patient.

6. During the operation, properly align the Ear Rods to the patient's ears so his/her head does not move.
7. Carefully tilt the patient's head back and adjust so his/her Frankfurt Line is vertical with the floor.
8. Direct the patient to swallow first before closing his/her mouth and to remain in his/her current position until image acquisition is completed.

9. Click the **READY** button on Console Software. No X-ray will be emitted at this point..
10. Now go to **4.3 X-ray Exposure** to start the exposure.

4.2.4 Waters' view Mode

Patient Positioning

1. Turn the Nasal Positioner to the **PA / Waters' view / Carpus** mode Positioning Marker as shown below.


2. Fold the Nasal Positioner up. The Nasal Positioner is not used in Waters' view mode.
3. Guide the patient to the CEPH unit.
4. Ask the patient to stand upright facing the X-ray detector. Make sure that the patient's shoulders are level and that his/her neck is relaxed.
5. Use the **Column UP/DOWN** button or switch option to adjust the height of the CEPH unit to approximately match the height of the patient.

WARNING

After adjusting the height of the column, align the Ear Rods to the patient.

6. During the operation, properly align the Ear Rods to the patient's ears so his/her head does not move.
7. Direct the patient to swallow first before closing his/her mouth, and guide the patient to bend the head backward 30° - 40°. Direct the patient to remain in the current position until image acquisition is completed.

8. Click the **READY** button on Console Software. No X-ray will be emitted at this point.
9. Now go to **4.3 X-ray Exposure** to start the exposure.

4.2.5 Carpus Mode


For Carpus Mode, install the Carpus Plate first before positioning the patient.

Installing the Carpus Plate

1. Turn the Nasal Positioner to the **PA / Waters' view / Carpus** mode Positioning Marker as shown below.


4. Fold the Nasal Positioner up. The Nasal Positioner is not used in Carpus mode.
5. Fit the two ends of the Carpus Plate into the two holes of the CEPH unit as below.

6. Confirm that the Carpus Plate is safely mounted.

Patient Positioning

1. Let the patient put his/her right hand splayed on the Carpus Plate as shown below. Make sure that the patient does not bend his/her fingers.


2. Ask the patient to close his/her eyes and stand still until the image acquisition is completed.
3. Click the **READY** button on Console Software. No X-ray will be emitted at this point.
4. Now go to **4.3 X-ray Exposure** to start the exposure.

4.3 X-ray Exposure

 WARNING	▪ If an emergency occurs during image acquisition, release the Exposure Switch to cease X-ray emission. ▪ The operator shall observe the X-ray safety regulations applicable to his/her area at all times during the operation of this equipment.
 CAUTION	▪ The operator must keep vocal/visual contact with the patient at all times during image acquisition process. ▪ Do not operate the PC during exposure. Doing so may cause the system to malfunction.
IMPORTANT	▪ Let the patient close the eyes during the operation. ▪ To acquire optimized images, instruct the patient to hold his/her breath and not to swallow.

1. Get out of the X-ray room and close the door.
2. Press and hold down the **Exposure Switch** until image acquisition is completed.

Yellow
- In Operation (X-ray On)

NOTICE	The image appears on the screen.
NOTICE	During X-ray exposure, the status appears as follows. ▪ The X-ray Exposure Indicator of the Exposure Switch and the LED light on the top of the equipment turn yellow. ▪ An alert sound comes out to indicate that X-ray emission is currently underway. ▪ On Console Software, the radiation mark turns yellow and "X-RAY" changes to "X-RAY ON".  X-RAY ON

3. Release the **Exposure Switch** when "Image capturing is completed" message appears on the screen.

4.4 Finishing the Scan

1. Leave enough space between the Ear Rods.
2. Fold the Nasal Positioner up in case it's unfolded.
3. Guide the patient out of the equipment.

4.5 Checking the Captured Images

Acquired images can be reconstructed and converted to DICOM format.

The exported images can be confirmed in **EzDent-i**.

NOTICE

Refer to the **EzDent-i** User Manual for more information.

1. The images are transferred to **EzDent-i** automatically.
2. The images are automatically saved if automatic save option is configured as default. If it is not configured as default, click **Save** button to save the images.
3. To check the image, double-click the one on the **Patient List**.

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The screenshot shows the EzDent-i main menu with the **PATIENT** tab selected. The interface includes a search bar, a recent patients list, and a table of patient records.

Chart No.	Name	Date of Birth
202191029_103203	Koo HyunJun	1974-08-15
20191029_103621	HyunJun Koo	1974-08-15
20191029_103722	Hyun koo	2019-01-01
202191029_103203	Koo HyunJun	1974-08-15

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The screenshot shows the **ADD PATIENT** dialog box. It contains the following fields:

- Chart No.*: 202191029_103203
- Name*: Last Name (Koo), First Name (HyunJun)
- Gender: Male
- Date of Birth: Year (1974), Month (8), Day (15)
- Doctor 1: None, Doctor 2: None
- Social ID: (empty)
- Phone: (empty)
- Mobile: (empty)
- E-Mail: (empty)
- Zip Code: (empty)
- Address: (empty)

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Main Menu
EzDent-i

PATIENT | ACQUISITION | VIEWER | CONSULT | REPORT

SEARCH
[Text Field]
[Search] [+]

RECENT PATIENTS LIST
[Recently Acquired]
[Recently Viewed]

DOCTOR
[All]

Chart No. 202191029_103203
Name Koo HyunJun
Gender/Age Male/47Y6M
Date of Birth 1974-08-15
[More Details]

Chart No.	Name	Date of Birth
202191029_103203	Koo HyunJun	1974-08-15
202191029_103621	HyunJun Koo	1974-08-15
202191029_103722	Hyun koo	2019-01-01
202191029_103203	Koo HyunJun	1974-08-15

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Date of Birth Year Month Day
1974 8 15

Doctor 1 None Doctor 2 None

Social ID
Phone
Mobile
E-Mail
Zip Code
Address

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
The screenshot shows the EzDent-i software interface. The top menu bar includes 'Main Menu', 'PATIENT', 'ACQUISITION', 'VIEWER', 'CONSULT', and 'REPORT'. The 'PATIENT' tab is active. On the left, there is a 'SEARCH' field with a 'Search' button and a '+' icon. Below it is a 'RECENT PATIENTS LIST' with a 'Recently Acquired' button. The main area displays patient information for 'Koo HyunJun' with a chart number of '202191029_103203', gender/age of 'Male/47Y6M', and date of birth of '1974-08-15'. A 'Date' button is visible on the right.

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The screenshot shows the EzDent-i software interface. The top menu bar includes 'Main Menu', 'PATIENT', 'ACQUISITION', 'VIEWER', 'CONSULT', and 'REPORT'. The 'PATIENT' tab is active. On the left, there is a 'SEARCH' field with a 'Search' button and a '+' icon. Below it is a 'RECENT PATIENTS LIST' with 'Recently Acquired' and 'Recently Viewed' buttons. The main area displays patient information for 'Koo HyunJun' with a chart number of '202191029_103203', gender/age of 'Male/47Y6M', and date of birth of '1974-08-15'. A 'More Details' button is visible. Below this, a table lists search results:

Chart No.	Name	Date of Birth
20191029_103203	Koo HyunJun	1974-08-15
20191029_103621	HyunJun Koo	1974-08-15
20191029_103722	Hyun koo	2019-01-01
202191029_103203	Koo HyunJun	1974-08-15

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The screenshot shows the EzDent-i main menu with the 'PATIENT' tab selected. On the left, there is a 'Main Menu' sidebar with options like 'SEARCH', 'RECENT PATIENTS LIST', and 'DOCTOR'. The main area displays patient information for 'Koo HyunJun' (Chart No. 202191029_103203, Date of Birth 1974-08-15). Below this, a table lists several patients, with 'HyunJun Koo' (Chart No. 20191029_103621) highlighted in green.

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
The screenshot shows the EzDent-i software interface. The top menu bar includes 'Main Menu', 'PATIENT', 'ACQUISITION', 'VIEWER', 'CONSULT', and 'REPORT'. The 'PATIENT' tab is active. On the left, there is a 'SEARCH' box with a search button and a '+' icon, and a 'RECENT PATIENTS LIST' with a 'Recently Acquired' button. The main area displays patient information for 'Koo HyunJun' with Chart No. 202191029_103203, Gender/Age Male/47Y6M, and Date of Birth 1974-08-15. A 'Date' button is visible on the right.

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The screenshot shows the EzDent-i software interface. The top menu bar includes 'Main Menu', 'PATIENT', 'ACQUISITION', 'VIEWER', 'CONSULT', and 'REPORT'. The 'PATIENT' tab is active. On the left, there is a 'SEARCH' box with a search button and a '+' icon, and a 'RECENT PATIENTS LIST' with 'Recently Acquired' and 'Recently Viewed' buttons. The main area displays patient information for 'Koo HyunJun' with Chart No. 202191029_103203, Gender/Age Male/47Y6M, and Date of Birth 1974-08-15. A 'More Details' button is visible. Below this, a table lists patient records:

Chart No.	Name	Date of Birth
20191029_103203	Koo HyunJun	1974-08-15
20191029_103621	HyunJun Koo	1974-08-15
20191029_103722	Hyun koo	2019-01-01
202191029_103203	Koo HyunJun	1974-08-15

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The screenshot shows the **ADD PATIENT** dialog box. The form includes fields for Chart No., Name, Gender, Date of Birth, Doctor 1, Doctor 2, Social ID, Phone, Mobile, E-Mail, Zip Code, and Address.

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
The screenshot shows the EzDent-i software interface. On the left is a 'Main Menu' with a 'SEARCH' section containing a text input field and a 'Search' button. Below it is a 'RECENT PATIENTS LIST' with 'Recently Acquired' and 'Recently Viewed' buttons. On the right, the search results for patient 'Koo HyunJun' are displayed, including 'Chart No. 202191029_103203', 'Name Koo HyunJun', 'Gender/Age Male/47Y6M', and 'Date of Birth 1974-08-15'. A 'More Details' button is visible below the patient information.

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The screenshot shows the EzDent-i software interface with search results. The 'SEARCH' section on the left is active. The results table below shows a list of patients with their 'Chart No.', 'Name', and 'Date of Birth'. The second row, '20191029_103621 HyunJun Koo 1974-08-15', is highlighted with a green border.

Chart No.	Name	Date of Birth
20191029_103203	Koo HyunJun	1974-08-15
20191029_103621	HyunJun Koo	1974-08-15
20191029_103722	Hyun koo	2019-01-01
202191029_103203	Koo HyunJun	1974-08-15

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The screenshot shows the EzDent-i main window with the 'PATIENT' tab selected. On the left is a 'Main Menu' with a search bar and buttons for 'Recently Acquired' and 'Recently Viewed'. Below that is a 'DOCTOR' dropdown set to 'All'. The main area displays patient information for 'Koo HyunJun' (Chart No. 202191029_103203, Date of Birth 1974-08-15). Below this is a table of patients:

Chart No.	Name	Date of Birth
202191029_103203	Koo HyunJun	1974-08-15
202191029_103621	HyunJun Koo	1974-08-15
202191029_103722	Hyun koo	2019-01-01
202191029_103203	Koo HyunJun	1974-08-15

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The 'ADD PATIENT' dialog box contains the following fields and controls:

- Chart No.*: 202191029_103203
- Name*: Last Name (Koo), First Name (HyunJun)
- Gender: Male (dropdown)
- Date of Birth: Year (1974), Month (8), Day (15)
- Doctor 1: None (dropdown), Doctor 2: None (dropdown)
- Social ID: [empty field]
- Phone: [empty field]
- Mobile: [empty field]
- E-Mail: [empty field]
- Zip Code: [empty field]
- Address: [empty field]
- Buttons: Add, Cancel

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
The screenshot shows the EzDent-i software interface. The top navigation bar includes tabs for PATIENT, ACQUISITION, VIEWER, CONSULT, and REPORT. The left sidebar has a 'Main Menu' with a search bar and buttons for 'RECENT PATIENTS LIST' and 'Recently Acquired'. The main area displays patient details for 'Koo HyunJun' with Chart No. 202191029_103203, Gender/Age Male/47Y6M, and Date of Birth 1974-08-15. A 'Date' filter is visible on the right.

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The screenshot shows the EzDent-i software interface with the search results list. The list contains four entries with their Chart No., Name, and Date of Birth. The entry with Chart No. 20191029_103621 and Name 'HyunJun Koo' is highlighted with a green border.

Chart No.	Name	Date of Birth
20191029_103203	Koo HyunJun	1974-08-15
20191029_103621	HyunJun Koo	1974-08-15
20191029_103722	Hyun koo	2019-01-01
202191029_103203	Koo HyunJun	1974-08-15

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Copyright © 2018 VATECH Co., Ltd.
Copyright by © 2018 VATECH Co., Ltd.